


E·Luminexx™

Vascular Stent



Clinical trial outcomes

A review of clinical trial data¹⁻⁴

The displayed clinical data are taken from a number of different studies, not a head-to-head study. Each study may have varying patient profiles and protocol structures that may affect the outcome rates set forth below, and therefore is not intended to demonstrate superiority of any one product.

The Luminexx™ Iliac Stent Clinical Study met the success criteria for the primary endpoint: posterior probability of at least 96% that the 9-month MACE rate was < 25%. MACE included peri-procedural death, TLR and stent segment restenosis, and non-study stent. (Non-study stent was added post hoc to the protocol definition of MACE.)

E·Luminexx™ Vascular Stent

(P080007)

94.03%

Primary
patency rate

98.72%

Anatomic success
(Site reposted)

0.75%

Mortality Rate
Devise or Procedure

The Luminexx™ Iliac Stent Clinical Study demonstrated competitive results when considering Primary Patency Rate, Anatomic Success and Mortality Rate.



Ordering information

BF Catheter System: 0.035" Guidewire				BF Catheter System: 0.035" Guidewire			
Order Code	Stent Diameter (mm)	Catheter Length (cm)	Stent Length (atm)	Order Code	Stent Diameter (mm)	Catheter Length (cm)	Stent Length (atm)
ZVM06020			20	ZVL09020			20
ZVM06030			30	ZVL09030			30
ZVM06040			40	ZVL09040			40
ZVM06050			50	ZVL09050			50
ZVM06060		80	60	ZVL09060	9	135	60
ZVM06080			80	ZVL09080			80
ZVM06100			100	ZVL09100			100
ZVM06120	6		120	ZVL09120			120
ZVL06020			20	ZVM10020			20
ZVL06030			30	ZVM10030			30
ZVL06040			40	ZVM10040			40
ZVL06050		135	50	ZVM10050		80	50
ZVL06060			60	ZVM10060			60
ZVL06080			80	ZVM10080			80
ZVL06100			100	ZVM10100			100
ZVL06120			120	ZVM10120			120
ZVM07020			20	ZVL10020	10		20
ZVM07030			30	ZVL10030			30
ZVM07040			40	ZVL10040			40
ZVM07050		80	50	ZVL10050		135	50
ZVM07060			60	ZVL10060			60
ZVM07080			80	ZVL10080			80
ZVM07100			100	ZVL10100			100
ZVM07120	7		120	ZVL10120			120
ZVL07020			20	ZVM12020			20
ZVL07030			30	ZVM12030			30
ZVL07040			40	ZVM12040			40
ZVL07050		135	50	ZVM12050		80	50
ZVL07060			60	ZVM12060			60
ZVL07080			80	ZVM12080			80
ZVL07100			100	ZVM12100			100
ZVL07120			120	ZVM12120	12		120
ZVM08020			20	ZVL12020			20
ZVM08030			30	ZVL12030			30
ZVM08040			40	ZVL12040			40
ZVM08050		80	50	ZVL12050		135	50
ZVM08060			60	ZVL12060			60
ZVM08080			80	ZVL12080			80
ZVM08100			100	ZVL12100			100
ZVM08120			120	ZVL12120			120
ZVL08020	8		20	ZVM14020			20
ZVL08030			30	ZVM14030			30
ZVL08040			40	ZVM14040			40
ZVL08050		135	50	ZVM14050		80	50
ZVL08060			60	ZVM14060			60
ZVL08080			80	ZVM14080			80
ZVL08100			100	ZVM14100			100
ZVL08120			120	ZVM14120	14		120
ZVM09020			20	ZVL14020			20
ZVM09030			30	ZVL14030			30
ZVM09040			40	ZVL14040			40
ZVM09050		80	50	ZVL14050		135	50
ZVM09060	9		60	ZVL14060			60
ZVM09080			80	ZVL14080			80
ZVM09100			100	ZVL14100			100
ZVM09120			120	ZVL14120			120

References

1. S.M.A.R.T.™ Control™ Nitinol Stent System. Premarket Approval Application Number P020036. Summary of Safety and Effectiveness Data, 2003. 2. Zilver® 635™ Vascular Stent. Premarket Approval Application Number P050017. Summary of Safety and Effectiveness Data, 2006. 3. Wallstent-Uni™ Iliac Endoprosthesis. Premarket Approval Application Number P940019. Summary of Safety and Effectiveness Data, 1996. 4. Ponc D, Jaff MR, Swischuk J, Feiring A, Laird J, Mehra M, Popma JJ, Donohoe D, Firth B, Heim E, Snead D. The nitinol SMART stent vs. Wallstent for suboptimal iliac artery angioplasty: CRISP-US trial results. J Vase Interv Radiol. 2004;15:911-918.

Please consult product labels and package inserts for indications, contraindications, hazards, warnings, cautions, and information for use



BD Switzerland

Sarl, Route de Crassier 17, Business Park Terre Bonne, Batiment A4, 1262 Eysins Switzerland,
Tel: +41 21 556 30 00. Fax: +41 44 722 5370

bd.com

