


EnCor™ MRI Probes and Introducer Set Specifications



Sample Notch Length	Measurement
Full sample	19 mm
Half sample	9.5 mm

Please consult product labels and inserts for any indications, contraindications, hazards, warnings, precautions and directions for use.

BD Switzerland Sarl
 Terre Bonne Park – A4, Route De Crassier, 17, 1262 Eysins, Vaud, Switzerland
 T: +41 21 556 3000

bd.com

BD, the BD Logo, EnCor and VisiLoc are trademarks of Becton, Dickinson and Company or its affiliates.
 © 2021 BD. All rights reserved. 0621/6168 BD-29892

CE
 2797

