


Dorado™

PTA Dilatation Catheter

Resistance
is futile



Confidence and precision for challenging cases

Enhanced Deliverability

- Checker™ Flex Points provide increased flexibility in tortuous anatomy.
- Improved shaft design promotes faster inflation and deflation.



Checker™ Flex Points

Longer lengths

- For treatment of long stenoses or post-dilatation of long stents with fewer inflations.

Balloon lengths up to 20cm




The Dorado™ PTA Dilatation Catheter provides strength and versatility with enhanced deliverability.


Resistance is futile

Ultra non-compliance

- Concentrates maximum dilatation force at the resistant lesion.
- Allows for inflation without the risk of overexpansion.
- Strong balloon material promotes consistent post-stent dilatation.
- Provides treatment of in-stent restenosis.



Semi-compliant balloon in a resistant stenosis.



Ultra non-compliant balloon in a resistant stenosis.

The Ultra non-compliant balloon is the optimal choice for challenging lesions and post-stent dilatation.

Ordering information

40 cm Shaft Length					
Order code	Balloon Diameter (mm)	Balloon Length (mm)	Nominal Pressure (atm)	Rated Burst Pressure (atm)	Sheath Size (F)
DR4042	4	20	8	24	5
DR4044		40	8	24	5
DR4052	5	20	8	24	5
DR4054		40	8	24	5
DR4062	6	20	8	24	6
DR4064		40	8	24	6
DR4072	7	20	8	22	6
DR4074		40	8	22	6
DR4082	8	20	8	22	6
DR4084		40	8	22	6
DR4094	9	40	8	20	6

80 cm Shaft Length					
Order code	Balloon Diameter (mm)	Balloon Length (mm)	Nominal Pressure (atm)	Rated Burst Pressure (atm)	Sheath Size (F)
DR8032	3	20	8	24	5
DR8034		40	8	24	5
DR80310		100	8	24	5
DR8042	4	20	8	24	5
DR8044		40	8	24	5
DR80410		100	8	24	5
DR8052	5	20	8	24	5
DR8054		40	8	24	5
DR8056		60	8	24	5
DR8058		80	8	24	5
DR80510		100	8	24	5
DR8062	6	20	8	24	6
DR8064		40	8	24	6
DR8066		60	8	24	6
DR8068		80	8	24	6
DR80610		100	8	22	6
DR80612		120	8	24	6
DR80615		150	8	24	6
DR80617		170	8	24	6
DR80620	200	8	24	6	
DR8072	7	20	8	22	6
DR8074		40	8	22	6
DR8076		60	8	22	6
DR8078		80	8	22	6
DR80710		100	8	22	6
DR80712		120	8	22	6
DR80715		150	8	22	6
DR80717	170	8	22	6	
DR80720	200	8	22	6	
DR8082	8	20	8	22	6
DR8084		40	8	22	6
DR8086		60	8	22	6
DR8088		80	8	22	6
DR80810		100	8	20	6
DR8092		9	20	8	20
DR8094	40		8	20	6
DR8098	80		8	20	6
DR80102	20		8	20	6
DR80104	10	40	8	20	6
DR80108		80	8	20	7

0.035" guidewire compatible

Units per case: 1

120 cm Shaft Length					
Order code	Balloon Diameter (mm)	Balloon Length (mm)	Nominal Pressure (atm)	Rated Burst Pressure (atm)	Sheath Size (F)
DR12044	4	40	8	24	5
DR120410		100	8	24	5
DR12052		20	8	24	5
DR12054	5	40	8	24	5
DR120510		100	8	24	5
DR12062	6	20	8	24	6
DR12064		40	8	24	6
DR120610		100	8	22	6
DR12074	7	40	8	22	6
DR120710		100	8	22	6
DR12084	8	40	8	22	6
DR120104	10	40	8	20	6

135 cm Shaft Length					
Order code	Balloon Diameter (mm)	Balloon Length (mm)	Nominal Pressure (atm)	Rated Burst Pressure (atm)	Sheath Size (F)
DR13532	3	20	8	24	5
DR13534		40	8	24	5
DR135310		100	8	24	5
DR13542	4	20	8	24	5
DR13544		40	8	24	5
DR135410		100	8	24	5
DR135412	4	120	8	24	6
DR135415		150	8	24	6
DR135417	5	170	8	24	6
DR135420		200	8	24	6
DR13552		20	8	24	5
DR13554	5	40	8	24	5
DR13556		60	8	24	5
DR13558		80	8	24	5
DR135510		100	8	24	5
DR135512		120	8	24	6
DR135515		150	8	24	6
DR135517		170	8	24	6
DR135520		200	8	24	6
DR13562	6	20	8	24	6
DR13564		40	8	24	6
DR13566		60	8	24	6
DR13568		80	8	24	6
DR135610		100	8	22	6
DR135612		120	8	24	6
DR135615	7	150	8	24	6
DR135617		170	8	24	6
DR135620		200	8	24	6
DR13572	7	20	8	22	6
DR13574		40	8	22	6
DR135710		100	8	22	6
DR135712		120	8	22	6
DR135715	8	150	8	22	6
DR135717		170	8	22	6
DR135720	8	200	8	22	6
DR13582		20	8	22	6
DR13584		40	8	22	6
DR135810	9	100	8	20	6
DR13594		40	8	20	6
DR135102	10	20	8	20	6
DR135104		40	8	20	6

0.035" guidewire compatible

Units per case: 1

Representatives name	Contact phone no.	Physicians signature
----------------------	-------------------	----------------------

BD Switzerland Sarl, Terre Bonne Park – A4, Route De Crassier, 17, 1262 Eysins, Vaud, Switzerland.

bd.com



Please consult product labels and instructions for indications, contraindications, hazards, warnings, and precautions.

BD, the BD logo, Dorado™ and all other trademarks are the property of Becton, Dickinson and Company or its affiliates. © 2021 BD. All rights reserved. Bard has joined BD. Illustrations by Mike Austin. Copyright © 2021. All rights reserved. BD-36485