

Frequently asked questions

1. Will I feel the suction used with the PureWick™ Urine Collection System?

No, the PureWick™ Urine Collection System is designed so that suction should not be felt. The device is held in place by your anatomy.

2. How often should the PureWick™ Female External Catheter be changed?

The device should be changed at least every 8–12 hours, or if soiled with faeces or blood.

3. Can the PureWick™ Urine Collection System be used when I'm on my side or sitting in a chair?

Yes – the PureWick™ Urine Collection System may be used when lying down or reclining or on your side, as well as when seated in a chair. Care must be taken to ensure that the product is not positioned beneath the tailbone.



Visit our website to learn more about the PureWick™ Urine Collection System or call **08000 90 24 93**

Script-easy® is the free continence and stoma home delivery service provided by Bard Limited. A Dispensing Appliance Contractor registered with the NHS, Script-easy® can dispense your continence and stoma products on receipt of your prescription - conveniently and discreetly delivered directly to your door whenever you need them. **Please note that the PureWick™ System is not available on prescription but Script-easy can still deliver it directly to your home free of charge.**

Indications for use: The PureWick™ Urine Collection System is to be used with the PureWick™ Female External Catheters which are intended for non-invasive urine output management.

Contraindications: The PureWick™ Female External Catheter should not be used on patients with urinary retention.

Please consult package insert for more detailed safety information and instructions for use.



Freephone:
08000 90 24 93

Email:
script-easy@bd.com

Nurse Advice Line Freephone:
0808 168 4048

Web:
www.scripteasy.co.uk

Becton, Dickinson U.K. Limited, 1030 Eskdale Road, Winnersh, Wokingham RG41 5TS, United Kingdom

www.purewickathome.co.uk

BD, the BD logo, PureWick and Script-easy are trademarks of Becton, Dickinson and Company. Script-easy® is the prescription dispensing service operated by Bard Limited. Bard is now a wholly-owned subsidiary of BD (Becton, Dickinson and Company). © 2022 BD and its subsidiaries. All rights reserved. BD-57066



PureWick™ Urine Collection System

Now available for home use.

A breakthrough for women with urinary incontinence.



What is the PureWick™ System?

The PureWick™ Urine Collection System is designed to better manage your urinary incontinence. It works with the PureWick™ Female External Catheter to wick urine away from the body into a designated collection canister. The PureWick™ Female External Catheter is changed every 8–12 hours and the collection canister containing urine is emptied into the toilet. You can now benefit from the PureWick™ Urine Collection System at home.

Benefits of using the PureWick™ System.



Helps keep the skin dry by wicking away urine.





May help to reduce the risk of nighttime falls, by staying in bed.*



Can make better sleep possible.*

How does the PureWick™ Urine Collection System work?



The PureWick™ Female External Catheter connects to the PureWick™ Urine Collection System.



The PureWick™ Female External Catheter is placed between your legs in order to collect urine.

***Based on survey responses (01/19) by 141 caregivers and users who purchased the PureWick™ System and a subsequent supply of PureWick™ Female External Catheters.