



Urological products

Product catalogue





Intermittent catheters	5
Magic3 Go™/Hydrosil Go™ Female.....	6
Magic3 Go™/Hydrosil Go™ Pediatric.....	7
Magic3 Go™/Hydrosil Go™ Male.....	8
Magic3™/Hydrosil™.....	9
Magic3™/Hydrosil™.....	10
BD™ Ready-To-Use.....	11
Closed systems and urethral trays	13
Magic3™ - Insertion Supply Kits	14
Incontinence solutions	17
PureWick™ Urine Collection System	18
PureWick™ Female External Catheter.....	19
Male external catheters	20
Spirit™ Hydrocolloid Adhesive Sheath	20
UltraFlex™ Self-Adhering Male External Catheter	21
Pop-On™ Self-Adhering Male External Catheter.....	21
WideBand™ Self-Adhering Male External Catheter	22
Natural™ Non-Adhesive, All-Silicone Male External Catheter.....	22
Self Adhesive silicon sheaths – Clear Advantage™.....	23
Transfix.....	23
Collection Systems	24
Uriplan™ Leg Bags (Sterile)	24
Uriplan™ Drainable Bed Bag (Sterile).....	24
Uriplan™ Non-Drainable Night Bag (Non-Sterile)	24
OneMT™ Single use Night Bags (Non-Sterile)	24
Spirit Leg Bags (Sterile).....	25
Spirit™ Drainable Bed Bag (Sterile).....	25
Spirit™ Single use Night Bag (Non-Sterile)	25
Prosys Leg Bag.....	25
Prosys Bed Bag	25
Urisac Leg Bag.....	26
Laprocare.....	26
Accessories	26
Essentials	27
Foley catheters and trays	29
Indwelling Foley Catheters.....	30
Bard™ Foley Trays	32
Bard™ Closed System Foley Trays.....	33







Intermittent catheters


Summary of indications and warnings

Intermittent catheters are intended for draining urine from the bladder. Reuse of a single use device may create a risk of patient or user infection. Please consult product labels and inserts for more detailed safety information and instructions for use.

Magic³ Go™ /Hydrosil Go™ Female

Hydrophilic Intermittent Catheter
with CleanGlide™ Technology

- Hydrophilic, female length catheter designed to glide easily and comfortably through the urinary tract
- Unique CleanGlide™ technology that maintains the catheter in a hydrated state so that it is ready to go out of the package
- Three-layered silicone catheter with an ultra-soft outer layer designed to help maximize comfort, firm middle layer for easy handling, and pliable innermost layer for easy navigation
- Comfort-designed tapered tip with four comfort-sized drainage eyes
- Robust handle for better control that is ergonomically designed for a secure grip that can easily be connected to a collection device
- Designed to minimize the risk of contamination by enabling the user to hold the catheter without touching the catheter shaft
- Discreet and compact consumer package that is easy to carry, easy to open, and easy to use right out of the package without the need of adding water or lubricant
- Environmentally conscious design - does not contain DEHP and/or phthalates



Female
Magic³ Go™
Female Intermittent Catheter

Product code	Length	French size	Qty
51808	6"	8 Fr.	30 /case
51810	6"	10 Fr.	
51812	6"	12 Fr.	
51814	6"	14 Fr.	
51816	6"	16 Fr.	
51818	6"	18 Fr.	

Female
Hydrosil Go™
Female Intermittent Catheter

Product code	Size Ch	Qty
71808	8	30 /case
71810	10	
71812	12	
71814	14	
71816	16	
71818	18	

Magic³ Go™ /Hydrosil Go™ Pediatric

Hydrophilic Intermittent Catheter
with CleanGlide™ Technology

100% silicone intermittent catheters

- Hydrophilic, pediatric length catheter designed to glide easily and comfortably through the urinary tract
- Unique CleanGlide™ technology that maintains the catheter in a hydrated state so that it is ready to go out of the package
- Three-layered silicone catheter with an ultra-soft outer layer designed to help maximize comfort, firm middle layer for easy handling, and pliable innermost layer for easy navigation
- Comfort-designed tapered tip with four comfort-sized drainage eyes
- Robust SureGrip™ insertion sleeve for better control that is ergonomically designed for a secure grip that can easily be connected to a collection device
- Designed to minimize the risk of contamination by enabling the user to hold the catheter without touching the catheter shaft
- Discreet and compact consumer package that is easy to carry, easy to open, and easy to use right out of the package without the need of adding water or lubricant
- Environmentally conscious design - does not contain DEHP and/or phthalates

Paediatric

Magic³ Go™

Paediatric Intermittent Catheter

Product code	Length	French size	Tip style	Qty/case
52806G	10"	6 Fr.	Straight	30/case
52808G	10"	8 Fr.	Straight	
52810G	10"	10 Fr.	Straight	

Paediatric

Hydrosil Go™

Paediatric Intermittent Catheter

Product code	Size Ch	Qty
72806G	6	30 /case
72808G	8	
72810G	10	



Magic³ Go™ /Hydrosil Go™

Male

Hydrophilic Intermittent Catheter
with CleanGlide™ Technology

- Hydrophilic, male length catheter designed to glide easily and comfortably through the urinary tract
- Unique CleanGlide™ technology that maintains the catheter in a hydrated state so that it is ready to go out of the package
- Three-layered silicone catheter with an ultra-soft outer layer designed to help maximize comfort, firm middle layer for easy handling, and pliable innermost layer for easy navigation
- Comfort-designed tapered tip with four comfort-sized drainage eyes.
- Robust SureGrip™ insertion sleeve for better control
- Designed to minimize the risk of contamination by enabling the user to hold the catheter without touching the catheter shaft
- Discreet and compact consumer package that is easy to carry, easy to open, and easy to use right out of the package without the need of adding water or lubricant
- Environmentally conscious design - does not contain DEHP and/or phthalates



Male
Magic³ Go™
Male Intermittent Catheter

Product code	Length	French size	Tip style	Qty
53810G	16"	10 Fr.	Straight	30/case
53812G	16"	12 Fr.	Straight	
53814G	16"	14 Fr.	Straight	
53816G	16"	16 Fr.	Straight	
53818G	16"	18 Fr.	Straight	
53820G	16"	20 Fr.	Straight	
50810G	16"	10 Fr.	Coudé	
50812G	16"	12 Fr.	Coudé	
50814G	16"	14 Fr.	Coudé	
50818G	16"	18 Fr.	Coudé	

Male
Hydrosil Go™
Male Intermittent Catheter

Product code	Size Ch	Qty
73810G	10	30 /case
73812G	12	
73814G	14	
73816G	16	
73818G	18	
73820G	20	

Magic³™/Hydrosil™

Hydrophilic Intermittent Catheter

- Ultra-soft outer layer designed to help maximize comfort, firm middle layer for easy handling, and pliable innermost layer for easy navigation
- Packaged with its own sterile water packet, immediate hydration creates virtually friction-free surface
- Comfort-designed tapered tip with four comfort-sized drainage eyes
- Environmentally conscious design - does not contain DEHP and/or phthalates

With SureGrip™ Insertion Sleeve

Product code	Length	French size	SureGrip™	Qty/case
52606G	10"	6 Fr.	Yes	30/Case
52608G	10"	8 Fr.	Yes	
52610G	10"	10 Fr.	Yes	
53610G	16"	10 Fr.	Yes	
53612G	16"	12 Fr.	Yes	
53614G	16"	14 Fr.	Yes	
53616G	16"	16 Fr.	Yes	
53618G	16"	18 Fr.	Yes	
53620G	16"	20 Fr.	Yes	



Product code	Length	French size	SureGrip™	Qty/case
51610	6"	10 Fr.	No	30/Case
51612	6"	12 Fr.	No	
51614	6"	14 Fr.	No	
51616	6"	16 Fr.	No	
51618	6"	18 Fr.	No	
52606	10"	6 Fr.	No	
52608	10"	8 Fr.	No	
52610	10"	10 Fr.	No	
53610	16"	10 Fr.	No	
53612	16"	12 Fr.	No	
53614	16"	14 Fr.	No	
53616	16"	16 Fr.	No	
53618	16"	18 Fr.	No	



Magic³™ /Hydrosil™

Silicone Intermittent Catheter

Intermittent Catheters	Product code	Size Ch	Qty
Coated Silicone Intermittent Catheters (Sterile)			
HydroSil® Gripper Male Hydrophilic Catheter	73610	10	30
	73612	12	30
	73614	14	30
	73616	16	30
	73618	18	30
HydroSil® Male Male Hydrophilic Catheter	63610M	10	30
	63612M	12	30
	63614M	14	30
	63616M	16	30
	63618M	18	30
	63620M	20	30
HydroSil® Rose Female Intermittent Catheter	71410	10	30
	71412	12	30
	71414	14	30
	71416	16	30
	71418	18	30

Intermittent Catheters	Product code	Size Ch	Qty
HydroSil® Female Female Hydrophilic Catheter	71610	10	30
	71612	12	30
	71614	14	30
	71616	16	30
	71618	18	30
HydroSil® Discreet Female Hydrophilic Catheter	61610C	10	30
	61612C	12	30
	61614C	14	30
	61616C	16	30
	61618C	18	30

Intermittent Catheters	Product code	Size Ch	Qty
HydroSil® Paediatric Gripper Hydrophilic Catheter	72606	6	30
	72608	8	30
	72610	10	30
HydroSil® Paediatric Hydrophilic Catheter	62606P	6	30
	62608P	8	30
	62610P	10	30

	Product code	Size Ch	Length	Qty
Coated PVC Intermittent Catheters (Sterile)				
Interglide® Coated Intermittent Catheter with finger grip (single use)	D6030	08–18	Male	25
	D6031	08–18	Female	25
	D6032	08–10	Paediatric	25
Uncoated PVC Intermittent Catheters (Sterile)				
Reliacath® with unique finger grip	D5030	12–16	Male	5
	D5031	12–16	Female	5
	D5032	08–10	Paediatric	5


BD™ Ready-To-Use

Hydrophilic Intermittent Catheter

- With an insertion aid for no-touch insertion technique including female and pediatric sizes
- PVP hydrophilic coating that is ready to use right out of the package
- Made with heat-polished eyelets designed for comfort
- Not made with DEHP or other phthalates
- Not made with BPA or natural rubber latex

Product code	French size	Length	Tip style	Units	HCPCS
RTU12M	12	Male - 16in	Straight	30	A4351
RTU14M	14	Male - 16in	Straight	30	A4351
RTU16M	16	Male - 16in	Straight	30	A4351
RTU06P	6	Pediatric - 10in	Straight	30	A4351
RTU08P	8	Pediatric - 10in	Straight	30	A4351
RTU10P	10	Pediatric - 10in	Straight	30	A4351
RTU10F	10	Female - 7.5in	Straight	30	A4351
RTU12F	12	Female - 7.5in	Straight	30	A4351
RTU14F	14	Female - 7.5in	Straight	30	A4351
Product code	French size	Length	Tip style	Units	HCPCS
RTU12C	12	Male - 16in	Coudé	30	A4352
RTU14C	14	Male - 16in	Coudé	30	A4352
RTU16C	16	Male - 16in	Coudé	30	A4352







Closed systems and urethral trays

Summary of indications and warnings

Intermittent catheters are intended for draining urine from the bladder. Reuse of a single use device may create a risk of patient or user infection. Please consult product labels and inserts for more detailed safety information and instructions for use.

Handle and dispose of in accordance with accepted medical practice and applicable local, state and federal laws and regulations.

Please consult product label and insert for any indications, contraindications, hazards, warnings, cautions and directions for use.

Magic³TM

Insertion Supply Kits

- With insertion supply kit, which includes Nitrile gloves, antiseptic wipes, underpad, and collection bag. 536XXGS codes include SureGripTM Insertion Sleeve

Product code	Length	French size	SureGrip TM	Qty/case
51614S	6"	14 Fr.	Yes	50/Case
53614GS	16"	14 Fr.	Yes	
53616GS	16"	16 Fr.	Yes	









Incontinence solutions

Summary of indications and warnings

The PureWick™ Urine Collection System is to be used with PureWick™ External Catheters which are intended for non-invasive urine output management. Do not use on individuals with urinary retention. Discontinue use if an allergic reaction occurs. Not recommended for users who are experiencing skin irritation or skin breakdown in device contact areas.

The PureWick™ Urine Collection System and accessories may be a potential biohazard. Clean and disinfect the PureWick™ Urine Collection System and accessories before disposal. Handle in accordance with accepted medical practice and applicable local and federal regulations.

Please consult the PureWick™ Urine Collection System user guide for more detailed safety information and instructions for use.

The PureWick™ Female External Catheters is intended for non-invasive urine output management in female patients. To avoid potential skin injury, never push or pull the PureWick™ Female External Catheter against the skin during placement or removal. Never insert into the vagina, anal canal, or other body cavities. Discontinue use if an allergic reaction occurs. Not recommended for patients who are: agitated, combative, or uncooperative and might remove the PureWick™ Female External Catheter.

After use, this product may be a potential biohazard. Dispose of in accordance with applicable local, state and federal laws and regulations.

Please consult the PureWick™ Female External Catheter insert for more detailed safety information and instructions for use.

Male external catheters are intended for use by male patients for the management of male urinary incontinence. Reuse of a single use device may create a risk of patient or user infection. Use of a single device for longer periods than 24 hours may increase the risk of complications. Do not use on irritated or compromised skin.

Handle and dispose of in accordance with accepted medical practice and applicable local, state and federal laws and regulations.

Please consult product label and insert for any indications, contraindications, hazards, warnings, cautions and directions for use.

PureWick™

Urine Collection System

- The PureWick™ Urine Collection System, the system to be used with the PureWick™ Female External Catheter, is an innovative option that is intended for non-invasive urine output management
- Designed to draw urine away, helping to keep skin dry¹



PureWick™ Urine Collection System

PureWick™ Female External Catheter

PureWick™ System

Product code	Description	Qty
PW200E	PureWick™ Urine Collection System	1
PWKIT03E	Accessory kit (<i>Replacement pump tubing, collector tubing and 2000cc collection canister</i>)	1

1. Ritter, J: Effectiveness and Functionality of a Female External Catheter. 2017

PureWick™

Female External Catheter

- Non-invasive way to manage urinary incontinence
- May reduce the risk of nighttime falls, by staying in bed¹
- Can make better sleep possible¹
- Can be placed in a few steps

Product code	Description	Qty
PWFX30E	PureWick™ Female External Catheter (<i>Not made with natural rubber latex</i>).	30

Attaches to the PureWick™
Urine Collection System or
existing hospital wall suction

Aligns with pubic
bone

Soft material wicks
urine
down to the
base

Continuous low pressure
suction at the base
removes urine



1. Based on survey responses (01/19) by 141 caregivers and users who purchased the PureWick™ System and a subsequent supply of PureWick™ Female External Catheters.



Spirit™

Hydrocolloid Adhesive Sheath

- Self-adhering male sheath to combine the excellent moisture-wicking properties of a hydrocolloid adhesive with the exceptional breathability of silicone¹
- A broad line of male external catheters with a diverse range of coverage
- Hydrocolloid adhesive wicks away moisture to help protect skin integrity, minimize skin maceration and maintain adhesion¹
- Made with 100% silicone




Style 1



Offers built-in adhesive in a length and placement found with traditional male external catheters.

Size		30/box	100/box	Size	30/box
Small	25mm	35301	35101	24mm	5243
Medium	29mm	35302	35102	28mm	5283
Intermediate	32mm	35303	35103	32mm	5323
Large	36mm	35304	35104	36mm	5363
X-Large	41mm	35305	35105	40mm	5403

Style 2



The 1.5" sheath with forward placement of adhesive designed to minimize unintended detachment.

Size		30/box	100/box	Size	30/box
Small	25mm	37301		24mm	7243
Medium	29mm	37302		28mm	7283
Intermediate	32mm	37303	37103	32mm	7323
Large	36mm	37304		36mm	7363
X-Large	41mm	37305		40mm	7403

Style 3




Our most popular style provides additional adhesive area compared to similar silicone self-adhering catheters. Combined with forward placement of adhesive designed to minimize unintended detachment.¹

Size		30/box	100/box	Size	30/box
Small	25mm	39301	39101	24mm	9243
Medium	29mm	39302	39302	28mm	9283
Intermediate	32mm	39303	39303	32mm	9323
Large	36mm	39304	39304	36mm	9363
X-Large	41mm	39305	39305	40mm	9403

1. Voegali, David, British Journal of Community Nursing, Vol. 18, No. 1, Jan. 2013

UltraFlex™

Self-Adhering Male External Catheter




- Offers built-in adhesive in a length and placement found with traditional male external catheters
- Made with 100% silicone

Size		30/box	100/box
Small	25mm	33301	33101
Medium	29mm	33302	33102
Intermediate	32mm	33303	33103
Large	36mm	33304	33104
X-Large	41mm	33305	33105



Pop-On™

Self-Adhering Male External Catheter



- The 1.5" sheath with forward placement of adhesive designed to minimize unintended detachment
- Made with 100% silicone

Size		30/box	100/box
Small	25mm	32301	32101
Medium	29mm	32302	32102
Intermediate	32mm	32303	32103
Large	36mm	32304	32104
X-Large	41mm	32305	32105



WideBand™

Self-Adhering Male External Catheter



- Provides additional adhesive area compared to similar silicone self-adhering catheters. Combined with forward placement of adhesive designed to minimize unintended detachment
- Made with 100% silicone

Size		30/box	100/box
Small	25mm	36301	36101
Medium	29mm	36302	36102
Intermediate	32mm	36303	36103
Large	36mm	36304	36104
X-Large	41mm	36305	36105



Natural™

Non-Adhesive, All-Silicone Male External Catheter



- Non-adhesive male external catheter. Each box comes with three reusable foam straps for securement
- Made with 100% silicone

Size		30 catheters and 3 straps/box
Small	25mm	38301
Medium	29mm	38302
Intermediate	32mm	38303
Large	36mm	38304
X-Large	41mm	38305



Male External Catheter

Self Adhesive silicon sheaths – Clear Advantage™

	Product code	Size mm	Qty
Style 1	786187	23	30
	786195	28	30
	786225	31	30
	786233	35	30
	786241	40	30

Transfix

	Product code	Size mm	Qty
Style 1	TF12530	25	30
	TF12930	29	30
	TF13230	32	30
	TF13630	36	30
	TF14130	41	30
Style 2	TF22530	25	30
	TF23930	29	30
	TF23230	32	30
	TF23630	36	30
	TF24130	41	30
Style 3	TF32530	25	30
	TF32930	29	30
	TF33230	32	30
	TF33630	36	30
	TF34130	41	30

***Caution:** This product contains natural rubber latex which may cause allergic reactions.



Collection Systems

- Made of heavy duty vinyl for long wear time
- Flip-Flo™ drainage valve opens and closes with a simple flip
- Reusable, sterile



Uriplan™ Leg Bags

Inlet Tube length	Product code	Vol	Qty
Direct inlet tube	BX3S/D3S	350ml	10
	BX5S/D5S	500ml	
	BX7S/D7S	750ml	
38cm inlet tube	BX7LX/D7LX	750ml	10
10cm inlet tube	BX5M	350ml	10
	D5M	500ml	
	BX7M/D7M	750ml	
30cm inlet tube	BX3L/D3L	350ml	10
	BX5L/D5L	500ml	
	BX7L/D7L	750ml	

Uriplan™ Drainable Bed Bag

Needle Free Sample Port | Non-Return Valve | 180° Tap

Inlet Tube length	Product code	Vol	Qty
98cm inlet tube	BX813131/ D813131	2000ml	10

Uriplan™ Non-Drainable Night Bag

Non-return Valve

Inlet Tube length	Product code	Vol	Qty
98cm inlet tube	D8420	2000ml	10

OneMT™ Single use Night Bags

Needle Free Sample Port | Non-Return Valve | 180° Tap

Inlet Tube length	Product code	Vol	Qty
90cm inlet tube	D1MT/D1MT30	2000ml	10

***Caution:** This product contains natural rubber latex which may cause allergic reactions.

Collection Systems

Spirit Leg Bags

Needle Free Sample Port | Non-return Valve | 180° Lever tap | 1 pair of Gloves per bag | 1 pair of Leg Straps

Inlet Tube length	Product code	Vol	Qty
10cm inlet tube	LB35010NS / LB35010S / DSLB350S	350ml	10
	LB50010S / DSLB500S	500ml	
	LB75010S / DSLB750S	750ml	
30cm inlet tube	LB35030S / LB35030NS / DSLB350L	350ml	10
	LB50030S / LB50030NS / DSLB500L	500ml	
	LB75030NS / LB75030S / DSLB750L	750ml	
60cm inlet tube	LB50060S / LB50060NS	500ml	10
	LB75060S / LB75060NS / DSLB750V	750ml	

Spirit™ Drainable Bed Bag

Needle Free Sample Port | Non-return Valve | 180° Lever tap | 1 pair of Gloves per bag

Inlet Tube length	Product code	Vol	Qty
120cm	BB1530NS	1500ml	10
120cm	BB2010NS	2000ml	10
120cm	BB2010S	2000ml	10
120cm	BB2030NS	2000ml	10
120cm	BB2030S	2000ml	10
120cm	DSBB2000	2000ml	10

Spirit™ Single use Night Bag T-tap

Inlet Tube length	Product code	Vol	Qty
90cm inlet tube	DSP1MT30	2000ml	30

Prosys Leg Bag

Inlet Tube length	Product code	Vol	Qty
Direct	PFORM500S / PFORM500V / PFORM500L / PFORM500VN	500ml	30
	PFORM750V	750ml	30
30cm inlet tube	PFORM500LN / PFORM500SN	500ml	30
	PFORM750SN / PFORM750VN	750ml	30

Prosys Bed Bag

Tube length	Product code	Vol	Qty
110cm	PFORM2000	2000ml	10
Bedhanger	PFORM2000B		1



Collection Systems

Urisac Leg Bag

Tube length	Product code	Vol	Qty
Direct	7660 / 7662 / 7663	350ml	

Laprocare

Tube length	Product code	Vol	Qty
90cm	LN103	1500ml	250
115cm	LN101 / LN102	2000ml	10
120cm	LN107	2000ml	20
120cm	LN108	4000ml	20

Accessories

Urisleeve™ Leg Bag Holder

Size	Product code	Length	Qty
Small	150111	24–39cm	4
Medium	150121	36–55cm	
Large	150131	40–70cm	
Extra Large	150141	65–95cm	

Uriplan™ Leg Bag Straps

Product code	Qty
15LS	5 pairs

PVC Bed Hanger for 2000ml bags

Product code	Qty
811602	1

Urilock™ Tap Lock

Product code	Qty
58081	10

Uristand™ Flat packed floor stand for Bed Bags

Product code	Qty
58081	10

***Caution:** This product contains natural rubber latex which may cause allergic reactions.

Essentials

Flip-Flo™ Catheter Valve

	Product code	Qty
Flip-Flo™ Catheter Valve	BFFS	5

Statlock™ Foley Stabilisation Device

	Product code	Qty
Statlock™ Foley Stabilisation Device	F0L0102DT	5

Catheter Retainer Straps

	Product code	Qty
Loc-Strap™ Catheter Retainer Strap	49999	5

Catheter Maintenance Solutions

	%	Product code	Length	Qty
Optiflo™ G	3.23% Citric Acid (Suby G)	CSG50	50ml	1
		CSG100	100ml	
Optiflo™ R	6.0% Citric Acid (Solution R)	CSR50	50ml	
		CSR100	100ml	
Optiflo™ S	0.9% Saline Solution	CSS50	50ml	
		CSS100	100ml	







Foley catheters and trays

Summary of indications and warnings

Foley catheters are intended for use in the drainage and/or collection and/or measurement of urine.

Please consult product label and insert for any indications, contraindications, hazards, warnings, cautions and directions for use.

Indwelling Foley Catheters

Long Term

Up to 12 weeks in situ, sterile

Tube length	Product code	Size Ch	Length	Balloon	Qty
Bard™ All Silicone	D165812-22	12-22	Standard	10ml	1
	D166112-16	12-16	Female	10ml	1

Hydrogel Coated All Silicone

Tube length	Product code	Size Ch	Length	Balloon	Qty
Lubri-Sil™ Aquafil™ Hydrogel Coated All Silicone Foley catheter (with pre-filled syringe of sterile water)	D175812-22E	12-22	Standard	10ml	1
	D176112-16E	12-16	Female	10ml	1
Lubri-Sil™ Aquafil™ Hydrogel Coated Foley catheter	D175808-10E	08-10	Paediatric	5ml	1

Hydrogel Coated Latex

Tube length	Product code	Size Ch	Length	Balloon	Qty
Biocath™ Aquamatic™ Hydrogel Coated Foley catheter (pre-filled with sterile water)	D226412-22	12-22	Standard	10ml	1
	D226812-22	12-22	Female	10ml	1
Biocath™ Hydrogel Coated Foley catheter	D226512-26	12-26	Standard	10ml	1
	D226612-26	12-26	Standard	30ml	1
	D226612-22	12-22	Female	10ml	1
	D226308-10	08-10	Paediatric	5ml	1



Indwelling Foley Catheters

Silicone Elastomer Coated Latex

Tube length	Product code	Size Ch	Length	Balloon	Qty
Silastic™ Foley catheter	33612-24CE	12–24	Standard	10ml	1
	33616-28CE	16–28	Standard	10ml	1

Short/Medium Term

Up to 28 days in situ, sterile

Tube length	Product code	Size Ch	Length	Balloon	Qty
BARDEX™ I.C. Acquafil™ Foley Catheter with Bactiguard® Silver Alloy Coating and Bard™ Hydrogel (with pre-filled syringe of sterile water)	D236512S-22S	12–22	Standard	10ml	1
	D236912S-16S	12–16	Female	10ml	1



Bardia™ Foley Trays



Long Term

Up to 12 weeks in situ, sterile

Biocath™ Bard™ Tray Preconnected to Uriplan™ Leg Bag

Inlet Tube	Product code	Size Ch	Length	Balloon	Qty
D5S Direct inlet	D22655S(12-16)	12-16	Standard	10ml	1
	D22695S(12-16)		Female		
D5M 10cm inlet	D22655M(12-16)	12-16	Standard	10ml	1
	D22695M(12-16)		Female		
D5L 30cm inlet	D22655L(12-16)	12-16	Standard	10ml	1
	D22695L(12-16)		Female		

Lubrisil™ Bard™ Tray Preconnected to Uriplan™ Leg Bag

Inlet Tube	Product code	Size Ch	Length	Balloon	Qty
D5S Direct inlet	D22655S(12-16)	12-16	Standard	10ml	1
	D22695S(12-16)	12-14	Female		
D5S 10cm inlet	D22655M(12-16)	12-16	Standard	10ml	1
	D22695M(12-16)	12-14	Female		
D5L 30cm inlet	D22655L(12-16)	12-16	Standard	10ml	1
	D22695L(12-16)	12-14	Female		

Short/Medium Term

Up to 28 days in situ, sterile

Bardex™ I.C. Bard™ Tray Preconnected to Uriplan™ Leg Bag

Inlet Tube	Product code	Size Ch	Length	Balloon	Qty
D5S Direct inlet	D23655S(12-16)	12-16	Standard	10ml	1
	D23695S(12-16)		Female		
D5M 10cm inlet	D23655M(12-16)	12-16	Standard	10ml	1
	D23695M(12-16)		Female		
D5L 30cm inlet	D23655L(12-16)	12-16	Standard	10ml	1
	D23695L(12-16)		Female		

Bardia™

Closed System Foley Trays

- 5cc silicone-elastomer coated latex Foley catheter pre-connected to a 2000ml urine drainage bag
- Waterproof underpad, fenestrated drape, 3 povidone iodine swabsticks, pre-filled 10cc inflation syringe, 5g lubricating jelly, 2 gloves not made with natural rubber latex, specimen container with label

Product code	French size	Qty/case
802016	16 Fr.	10/case
802018	18 Fr.	

Caution: This product contains natural rubber latex which may cause allergic reactions.









For more information on any of our products,
visit www.bd.com.

BD Switzerland Sarl, Terre-Bonne Park - A4,
Route de Crassier 17, 1262 Eysins, Switzerland

bd.com

BD, the BD Logo and Bard, Bardia, Care, Comfort, Choice., Clean-Cath, CleanGlide, Dispoz-a-Bag, Flip-Flo, Magic¹, Magic³ Go, Natural, Personal Catheter, Pop-On, PureWick Spirit, StatLock, Sure-Grip, Touchless, UltraFlex and WideBand are trademarks of Becton, Dickinson and Company and its affiliates.
© 2023 BD. All rights reserved. 0906-23 BD-92040

CE



BD